

The mini homogenizer MT-13K is hand-held designed. It is smart and light, easy to operate. With DC motor to drive the grinding pestle, configure 1.5ml centrifugal tube, efficiently grind trace tissue. The grinding pestle is light and reusable, suitable to a variety of sterile ways, convenient to be dismantled.

【 Applications 】

1. Suitable for grinding various plant tissue sample such as root, stem, leaf, flower, fruit, seed, etc.
2. Suitable for grinding various animal tissue sample such as cerebrum, heart, lung, stomach, liver, thymus, kidney, intestine, lymphonodus, muscle, skeleton, etc.
3. Suitable for homogenate.

【 Product features 】

1. Light and easy to operate, broad scope of application.
2. Use micro DC motor for grinding trace tissue, efficient and effective.
3. Grinding pestle is light and reusable, suitable for various sterile ways, convenient to be dismantled.

【 Technical Parameters 】

Type	MT-13K Mini Hand Held Homogenizer
Speed	13000rpm
Power	AC 220V
Noise	60dB
Net Weight	0.1kg
Dimension	25 x 25 x 160 mm

